

A-Tech Products

SQA-V for Animals

The A-Tech animal line of SQA's combine a fast, precise and accurate sperm quality analyzer with PC dosing and data management software. This winning combination allows the laboratory to test FRESH semen and complete the dosing process in about 3 minutes. Each SQA is species specific: BULL, EQUINE, PIG and TURKEY and reports values and presents data that is relevant to the industry.

Features:

- User friendly, menu driven testing and dosing instructions
- Multi use capillaries reduce cost per test
- Standardizes and enhances the quality and through-put of the testing and dosing process
- Includes on-board printer and microscope (x300 & x500)



SQA-Vb BULL

Tests FRESH and FROZEN samples:

- Total Sperm Concentration (TSC)
- Motile Sperm Concentration (MSC)
- Progressive Motile Sperm Concentration (PMSC)
- % Motility
- Morphology (FRESH ONLY)
- Velocity



SQA-Ve EQUINE

Tests RAW, EXTENDED, COOLED and FROZEN samples:

- Total Sperm Concentration (TSC)
- Motile Sperm Concentration (MSC)
- Progressive Motile Sperm Concentration (PMSC)
- % Motility
- % Progressive Motility
- Velocity
- % Morphology (RAW only)



SQA-Vp PIG

Tests both FRESH and EXTENDED samples:

- Total Sperm Concentration (TSC)
- Motile Sperm Concentration (MSC)
- Motility Grading
- Morphology (FRESH only)



SQA-Vt TURKEY

Tests NEAT and DILUTED pooled samples:

- Total Sperm Concentration (TSC)
- Motile Sperm Concentration (MSC)
- % Motility

Dosing and Data Management Software

A PC based dosing and data management software program is included with each SQA-V sold.

Features:

- Automatically dose FRESH/RAW samples tested on the SQA by following clear diluting/extending instructions
- View/Print reports, graphs and histograms
- Optional language settings
- Windows XP and VISTA compatible

www.a-tech-global.com

