

QwikCHECK™ GOLD PIG

Sperm Quality Analyzer

Measure both FRESH and EXTENDED BOAR semen on the farm in less than ONE minute!

Low cost, simple and accurate!



- Test FRESH/EXTENDED BOAR semen
- Print test results to a label maker
- On-board heating system for sample testing
- Sample requirement: 500 µl
- Re-usable, washable testing capillary

FRESH and EXTENDED SEMEN tests performed:

- Total Sperm Concentration M/ml
- % Motility
- Motility Grading # 1-4
- Motile Sperm Concentration
- Total # Sperm
- Total # Motile Sperm



“Remember, it all started with a Sperm!”